

# Sample Stand



***Use the same sensor in the laboratory, process scale-up, and manufacturing.***

This approach gives reproducible results and makes scale up from development to production easy. The analyzer is easy to move and remove from the sample stand, blender, dryer, or coater.

## **Specifications**

Materials: All stainless steel construction

Probe fitting: 3" sanitary fitting with tri-clover clamp

Probe fitting dimension: ID 2.75" round x 2"

Base dimension: 8" square x 3 1/2" H

**PART NUMBER: SS-100**

Control Development, Inc., reserves the right to change specifications without notice.

Control Development, Inc.  
2633 Foundation Drive, South Bend, IN 46628  
Phone: (574) 288-7338 Fax: (574) 288-7339  
[www.controldevelopment.com](http://www.controldevelopment.com)  
sales@controldevelopment.com

**CONTROL  
DEVELOPMENT**  
THE PRICE PERFORMANCE LEADER...